

”More independence
with your personal care”



ROPOX TOILET SUPPORT ARMS

GET SOMETHING TO LEAN ON

Toilet support arms are a great aid if you need to get on and off the toilet safely. Their ergonomic shapes provide a comfortable grip, and the cleaning-friendly surfaces optimize hygiene.

Because they are designed with two heights and an anti-slip surface, you have a firm grip from both standing and seated position. If you need a little extra support for getting up and sitting down, you can mount the standing/sitting attachments.

The support arms can be adjusted in height and at different angles to suit the individual user's needs.

If you need extra room in the bathroom, the support arms can be folded up to vertical position.






ROPOX™
BETTER WAYS TO BETTER DAYS
- A PART OF AddLife®

ROPOX TOILET SUPPORT ARMS




Height adjustment: 12 cm
 Max. user weight: Tested acc. to DS/ISO 17966:2016
 (see max. user weight below)

VARIANT/MODEL

WAVE - 90 CM		
		
190 kg Item no.: 40-40115	390 kg Item no.: 40-40125	385 kg Item no.: 40-40145

WAVE - 76 CM

255 kg Item no.: 40-40151


STRAIGHT - 90 CM		
		
190 kg Item no.: 40-40110	390 kg Item no.: 40-40120	385 kg Item no.: 40-40140

STRAIGHT - 76 CM


285 kg/385 kg Item no.: 40-40150/40-40152

ACCESSORIES

Height adjustment
 (adjusted with tool)
 Item no.: 40-40906



Height adjustment
 (adjusted without tool)
 Item no.: 40-40907



Toilet roll holder
 Item no.: 40-40910



Standing/Sitting attachment
 Item no.: 40-40925



NB! Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.

The items are produced according to the following Standard: DS / ISO 17966:2016

ROPOX A/S · Ringstedgade 221 · 4700 Naestved, Denmark – Phone: +45 5575 0500 · info@ropox.com · www.ropox.com

ROPOX™
 BETTER WAYS TO BETTER DAYS
 – A PART OF AddLife®

No. 1962390-1122